


AvanCE[®] capsule endoscope delivery device




Direct placement of video capsules...

The **AdvanCE capsule endoscope delivery device** allows for direct endoscopic placement of video capsules in patients who are unable to swallow or pass the capsule through the pylorus.

The AdvanCE device...

- Facilitates endoscopic placement of video capsules into the stomach or duodenum
- Minimizes risk of stalled capsule in esophagus or stomach, related to anatomical or motility disorders.



For more information on the AdvanCE delivery device, contact STERIS Endoscopy today.

product number	sheath diameter (mm)	length (cm)	compatible video capsule dimensions (mm)	units/box
BX00711144	2.5	180	diameter: 10.5 - 11.5 length: 23.5 - 26.5	3

Distributed by:

 **STERIS**
Endoscopy

STERIS
5960 Heisley Road
Mentor, OH 44060-1834 – USA
Customer Service: +1 800-548-4873
Fax: 440-639-4495
www.steris.com

Olympus® is a registered trademark of Olympus Corporation.
© 2020 STERIS Endoscopy or its affiliates. All rights reserved.
Unless otherwise indicated, all marks denoted with ® or ™ are registered with the U.S. Patent and Trademark Office,
or trademarks owned by STERIS or its affiliates.

Follow us:    

760398 Rev. E